




## Gastro-Intestinal Endoscopes

### Features

- Insertion Tube Diameter: 7.9 mm or 9.8 mm
- Working Length: 1.5 m, 2.0 m or 2.7 m
- Total Length: 1.75m, 2.25m or 3.0m
- Instrument Channel Diameter: 2.0 mm or 2.8 mm
- Air/Water, Suction and Instrument Channels
- 4-way Tip Articulation



**Package includes 150W Halogen Light Source, Water Tank and Air Pump, Biopsy Forceps, Cleaning Brushes, Leakage Tester, Jaw Spacers, Water Proof Insertion Tube and a Carrying Case**

### Medit Flexible GI-Endoscopes

The Medit Gastro-Intestinal Endoscopes are designed to perform diagnostic procedures and post operative examinations in the stomach and colon of cats, dogs and horses. These flexible Fiberscopes feature a narrow 7.9 mm insertion tube for use in cats and dogs or a larger, 9.8 mm insertion tube for large breed dogs and horses. Both models feature 20,000 fiber optic pixels, a waterproof, coated insertion tube with 4-way articulation and channels for air/water, suction and instruments. As with all of Medit Inc.'s products, we also offer video and image capture systems, which will allow you to record and save images and video for future reference.



For more information or to place your order contact:

### MEDIT INC. VETERINARY ENDOSCOPES

96 Sherbrook Street Winnipeg, MB, R3C 2B3, Canada [www.endoscopy4you.com](http://www.endoscopy4you.com)

Toll Free Phone: (800) 239-9934 Toll Free Fax: (888) 686-8600

Email: [office@meditinc.com](mailto:office@meditinc.com)

